## MAYFIELD® Adult Radiolucent Horsehoe Headrest ( **REF** - A2010)

### Instruction Manual

EN – ENGLISH .....	2
FR – FRANÇAIS .....	16
IT – ITALIANO .....	30
DE – DEUTSCH .....	44
ES – ESPAÑOL .....	58
NL – NEDERLANDS .....	72

Manufacturer:

Integra LifeSciences Corporation  
4900 Charlemar Drive, Building A  
Cincinnati, OH 45227, USA  
Tel: 513-533-7979  
Fax: 513-271-1915  
integralife.com

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# Rx ONLY

This page is intentionally left blank.

**MAYFIELD® Adult Radiolucent  
Horseshoe Headrest  
( REF - A2010)**

**Instruction Manual**

## Meaning Of Symbols Used In This Manual - ENGLISH

 CAUTION! 

Hazards which could result in equipment or property damage

 WARNING! 

Hazards which could result in severe personal injury or death

Caution

Product complies with the requirements of MDR 2017/745

Manufacturer

Authorized Representative in the European Community

Consult Instructions for Use

**Rx ONLY**

Caution: Federal (US) law restricts this device to sale by or on the order of a licensed healthcare practitioner.

Catalog number

Date of manufacture (YYYY-MM-DD)

Lot number

Serial number

Medical Device

## EN – ENGLISH

### Indications for Use/Intended Purpose

The MAYFIELD® Adult Radiolucent Horseshoe Headrest is intended to be used to support a patient's head in both the prone and supine positions during diagnostic examination or surgical procedures where support and stabilization of the head is desired but complete fixation of the head is not required.

**WARNING:**

Failure to read and follow instructions furnished in this product insert may result in serious patient injury..

**WARNING:**

Device must not be used if it appears to be damaged or functioning incorrectly. Doing so may result in serious patient injury.

**WARNING:**

Ensure the mating starburst teeth of the Swivel Adaptor and the Stationary Horseshoe Faceplate are the same size and have properly engaged (meshed) when fully tightened. Failure to do so may damage teeth, compromise stability, and result in serious patient injury.

**WARNING:**

Do not position the Slide Horseshoe Faceplate past the end of the Stationary Horseshoe Faceplate slide bar (Figure 11). Doing so may result in serious patient injury.

**WARNING:**

Horseshoe Faceplates should be adjusted prior to patient positioning. Make certain to always support the patient's head when adjusting the Horseshoe Faceplates or unlocking any connection of this device. Failure to do so may result in serious patient injury.

**WARNING:**

Failure to properly position the patient and to fully secure all adjustment portions of this or any support device may result in serious patient injury.

**WARNING:**

Pulley Rod may be used for traction of up to 40 lbs. Exceeding 40 lbs. of traction load may result in device failure and serious patient injury.

**WARNING:**

Ensure the mating starburst teeth of the Pulley Rod and the Stationary Horseshoe Faceplate are the same size and have properly engaged (meshed) when fully tightened. Failure to do so may damage teeth, compromise stability and support of traction load, and result in serious patient injury.

**WARNING:**

Device must not be used if it has been dropped or mishandled. Doing so may result in serious patient injury.

**WARNING:**

Do not steam sterilize any component of this device. Doing so may damage device and result in serious patient injury.

## Description

The MAYFIELD Adult Radiolucent Horseshoe Headrest (Figure 1) includes:

1. Stationary Horseshoe Faceplate
2. Slide Horseshoe Faceplate
3. Adjustment Knob
4. Pulley Rod
5. Pulley Rod Knob
6. Horseshoe Gel Pad, Adult Left
7. Horseshoe Gel Pad, Adult Right


Figure 1

The MAYFIELD Adult Radiolucent Horseshoe Headrest is designed to provide support and stabilization of the head in conjunction with the following equipment: A1079 MAYFIELD Radiolucent Base Unit (shown in Figure 2), A2004 MAYFIELD Radiolucent Headrest System, A1111 MAYFIELD Radiolucent Tri-Star Swivel Adaptor, A2079 MAYFIELD Infinity XR2 Base Unit, and A2111 MAYFIELD Infinity XR2 Tri-Star Swivel Adaptor.


Figure 2

The MAYFIELD Adult Radiolucent Horseshoe Headrest is designed to attach to the Short Swivel Adaptor (Figure 3) of the MAYFIELD Radiolucent Base Unit or the MAYFIELD Infinity XR2 Base Unit (not shown). If extra starburst attachments for ancillary equipment are desired, the MAYFIELD Adult Radiolucent Horseshoe Headrest may be attached to the MAYFIELD Radiolucent Tri-Star Swivel Adaptor (REF A-1111, Figure 4) of the MAYFIELD Radiolucent Headrest System (REF A-2004) in conjunction with the MAYFIELD Radiolucent Base Unit or the MAYFIELD XR2 Tri-Star Swivel Adaptor (REF A2111 (not shown)).


Figure 3


Figure 4

The MAYFIELD Adult Radiolucent Horseshoe Headrest is designed to be used with the included Horseshoe Gel Pads. Replacement Gel Pads are available individually (REF 41C1444 Adult Left and 41C1445 Adult Right).

The included Pulley Rod can be attached for traction of up to 40 lbs. if desired, or left unattached.

Use of MAYFIELD products and accessories in conjunction with other manufacturers' stabilization equipment is not recommended.

### **Intended Population**

MAYFIELD Skull Clamp fixation devices are not recommended for use on children under five (5) years of age. Extreme caution should be exercised in pediatric cases because of the thin skull.

### **Inspection**

Always inspect instruments before and after use. If a component appears damaged and/or does not seem to function properly, do not use the device and immediately send the instrument to an authorized Integra repair center for evaluation, repair or replacement. Allow your Integra Representative to inspect this device a minimum of two times per year to assist you in its proper function.

#### **WARNING:**

Device must not be used if it appears to be damaged or functioning incorrectly. Doing so may result in serious patient injury.

### **Gel Pads**

1. Pads should not have any rips or tears in the top surface.
2. Backs of the pads should not be broken or cracked.

### **Stationary Horseshoe Faceplate with attached Pulley Rod Knob**

1. Pulley Rod Knob should freely turn and easily screw into the threaded hole of the Pulley Rod.
2. Small starburst teeth of the Stationary Horseshoe Faceplate should match up and tighten securely with the starburst teeth of the Pulley Rod.
3. Large starburst teeth of the Stationary Horseshoe Faceplate should match up and tighten securely with the starburst teeth of the Radiolucent Short Swivel Adaptor or Radiolucent Tri-Star Swivel Adaptor.
4. Parts should be smooth and free of cracks and voids.
5. Gel Pad should slide easily onto the Stationary Horseshoe Faceplate and snap in place.

### **Slide Horseshoe Faceplate with attached Adjustment Knob**

1. Slide Horseshoe Faceplate should slide easily onto the Stationary Horseshoe Faceplate.
2. Adjustment Knob should turn easily and securely lock down on the slide bar of the Stationary Horseshoe Faceplate.
3. Parts should be smooth and free of cracks and voids.
4. Gel Pad should slide easily onto the Slide Horseshoe Faceplate and snap into place.

### **Pulley Rod**

1. Part should be smooth and free of cracks and voids.

## **Instructions for Use**

### **Mounting of Base Unit and attachment of Transitional Members and Swivel Adaptor**

#### **IMPORTANT:**

Refer to the corresponding MAYFIELD instruction manual (REF 451A2004, 451A2079) for the correct usage of the applicable MAYFIELD Base Unit prior to using this device. If applicable, refer to the corresponding MAYFIELD Tri-Star Swivel Adaptor instruction manual (REF 451A1111, 451A2111) prior to using this device.

1. Mount the MAYFIELD Radiolucent Base Unit to the operating table.
2. Attach the Transitional Members and the Short Swivel Adaptor to the Radiolucent Base Unit. If extra starburst attachments for ancillary equipment are desired, the MAYFIELD Radiolucent Tri-Star Swivel Adaptor may be used instead of the Short Swivel Adaptor.

### Attachment of Adult Radiolucent Horseshoe Headrest to Swivel Adaptor

1. Attach the Adult Radiolucent Horseshoe Headrest to the chosen Swivel Adaptor by inserting the torque screw of the Short Swivel Adaptor (Figure 5) or the torque screw of the Radiolucent Tri-Star Swivel Adaptor (Figure 6) into the threaded hole of the Stationary Horseshoe Faceplate such that the starburst teeth of the Swivel Adaptor are facing the large starburst teeth of the Stationary Horseshoe Faceplate.


Figure 5


Figure 6

2. Turn the Swivel Adaptor torque screw clockwise several times.
3. Position the Adult Radiolucent Horseshoe Headrest in desired orientation and turn the Swivel Adaptor torque screw clockwise until fully tightened (Figures 7 and 8)


Figure 7


Figure 8

#### WARNING:

Ensure the mating starburst teeth of the Swivel Adaptor and the Stationary Horseshoe Faceplate are the same size and have properly engaged (meshed) when fully tightened. Failure to do so may damage teeth, compromise stability, and result in serious patient injury.

4. Slide the Gel Pads onto the Horseshoe Faceplates until snapped in place.

### Adjustment of Horseshoe Faceplates

1. Turn the Adjustment Knob counterclockwise until the Slide Horseshoe Faceplate can slide freely (Figure 9).


Figure 9

2. Position the Slide Horseshoe Faceplate with respect to the Stationary Horseshoe Faceplate to accommodate the size of the patient's head.
3. Turn the Adjustment Knob clockwise until fully tightened (Figure 10).



Figure 10

**WARNING:**

Do not position the Slide Horseshoe Faceplate past the end of the Stationary Horseshoe Faceplate slide bar (Figure 11). Doing so may result in serious patient injury.

MAXIMUM POSITION  
WHEN FACES FLUSH

Figure 11

**WARNING:**

Horseshoe Faceplates should be adjusted prior to patient positioning. Make certain to always support the patient's head when adjusting the Horseshoe Faceplates or unlocking any connection of this device. Failure to do so may result in serious patient injury.

**WARNING:**

Failure to properly position the patient and secure all connections on this device may result in serious patient injury.

### Attachment of Pulley Rod

**WARNING:**

Pulley Rod may be used for traction of up to 40 lbs. Exceeding 40 lbs. of traction load may result in device failure and serious patient injury.

1. Insert threaded portion of the Pulley Rod Knob into the threaded hole of the Pulley Rod such that the starburst teeth of the Pulley Rod are facing the small starburst teeth of the Stationary Horseshoe Faceplate (Figure 12).


Figure 12

2. Turn the Pulley Rod Knob clockwise several turns.
3. Position the Pulley Rod in desired orientation.
4. Turn the Pulley Rod Knob clockwise until fully tightened (Figure 13).


Figure 13

**WARNING:**

Ensure the mating starburst teeth of the Pulley Rod and the Stationary Horseshoe Faceplate are the same size and have properly engaged (meshed) when fully tightened. Failure to do so may damage teeth, compromise stability and support of traction load, and result in serious patient injury.

## Cleaning

### Group 1 - Non-Critical<sup>1</sup>, Heat Sensitive, pH Sensitive

- Headrests

### Group 2 - Non-Critical<sup>1</sup>, Heat Sensitive, pH Sensitive, Alcohol Sensitive

- Gel Pads

**WARNING:**

There is no validated automated cleaning/washing/disinfection method for Group 2 (gel pads). A thorough, manual cleaning process is recommended for these devices. Automated cleaning methods may not be effective and may result in damage and reduced performance.

Prior to cleaning, disassemble components.

1. Remove the Gel Pads from the Horseshoe Faceplates.
2. Turn the Adjustment Knob counterclockwise until the Slide Horseshoe Faceplate can freely separate from the Stationary Horseshoe Faceplate.
3. Turn the Pulley Rod Knob counterclockwise until the Pulley Rod can freely separate from the Stationary Horseshoe Faceplate.

**CAUTION:**

The Adjustment Knob is permanently attached to the Slide Horseshoe Faceplate. Use of excessive force in an attempt to remove the Adjustment Knob from the Slide Horseshoe Faceplate can result in damage to components.

**CAUTION:**

The Pulley Rod Knob is permanently attached to the Stationary Horseshoe Faceplate. Use of excessive force in an attempt to remove the Pulley Rod Knob from the Stationary Horseshoe Faceplate can result in damage to components.

The MAYFIELD Adult Radiolucent Horseshoe Headrest should be thoroughly cleaned after each use. Clean all components, including the Gel Pads, with a neutral pH detergent. A soft brush can be used if necessary. Clean all components thoroughly to prevent traces of blood or debris from interfering with proper function of the device.

Rinse thoroughly with de-ionized or sterile water to remove any residual detergent.

All components of the MAYFIELD Adult Radiolucent Horseshoe Headrest except the Gel Pads may be wiped down with a 70% isopropyl alcohol solution.

**WARNING:**

Do not steam sterilize any component of this device. Doing so may damage device and result in serious patient injury.

**CAUTION:**

Do not use an alcohol solution on the Gel Pads, immerse the Gel Pads in liquid, or expose the Gel Pads to heat via steam sterilization or any other means.

## Manual Wash

### CAUTIONS

- Alkaline and highly acidic detergents and solutions cause damage to the devices.
- Channels and crevices found on this device require particular attention during cleaning.
- Pay special attention to the water quality used throughout reprocessing. Hard water can damage the surface of the equipment. Avoid using hard water. Instead use purified water unless otherwise specified.

### Limitations on reprocessing

- Repeated processing has minimal effects on these devices. Product life is normally determined by wear and damage due to use.
- It is important to have Integra NeuroSpecialists perform routine inspections (twice yearly is recommended). See contact information below.

## INSTRUCTIONS

### Containment/Transportation

- Follow health care facility protocol for safe containment and transport to the decontamination environment.
- It is recommended that devices are cleaned immediately after use.

---

<sup>1</sup> AAMI TIR 12:2010 Designing, testing, and labeling reusable medical devices for reprocessing in Health Care facilities.

**Preparation for Cleaning**

- Disassembly is required.

Disassembling Equipment	Inspecting for Cleanliness
Remove the gel pads from the unit by sliding off. Keep all joints in the unlocked position. Slide the adjustable headrest base until it comes apart.	Inspect the teeth of the starbursts and make sure to inspect all joints. Pay close attention to the area around the locking screw.

**Cleaning – Manual Equipment: Water, Neutral pH Detergent, Soft Bristle Brush, Towels**

**Method**

1. Prepare neutral pH enzymatic detergent solution (e.g. Endozime® AW Triple Plus with APA (Ruhof), 1:128 ratio) according to detergent manufacturer's instructions using lukewarm tap water.
2. Prepare equipment for soaking by disassembling removable parts and loosening connections.
3. Rinse equipment in warm water before placing into bath.
4. Completely soak equipment in water/detergent solution for 30 minutes maximum.
5. Clean thoroughly with a soft nylon bristle brush. NOTE: If possible, use a disposable brush.
6. Rinse in warm purified water until all visible substances and residual detergent are removed. NOTE: Make sure to give special attention to hard-to-reach areas.
7. Thoroughly dry equipment with soft clean towels and use medically compressed air if needed, to dry channels, crevices and lumens.
8. Inspect the equipment to make sure there is no visible organic debris or residue from cleaning agent.

Repeat process if any soil is detected.

**Drying**

- Products should be dry at this point. If wetness or excess liquid is detected, dry with a soft clean towel.
- Medically compressed air can be used if needed.

**Optional Automatic Wash / Disinfect**

**CAUTIONS**

- Alkaline and highly acidic detergents and solutions cause damage to the devices.
- Channels and crevices found on these devices require particular attention during cleaning.
- Pay special attention to the water quality used throughout reprocessing. Hard water can damage the surface of the equipment. Avoid using hard water. Instead use purified water unless otherwise specified.

**Limitations on Reprocessing**

- **Repeated processing has negative effects on these devices and is not recommended for routine use.**
- It is important to have Integra NeuroSpecialists perform routine inspections (twice yearly is recommended). See contact information.

**INSTRUCTIONS**

**Cleaning - Automated Equipment: Neutral pH Detergent**

**Method**

1. Prepare equipment for cleaning by disassembling removable parts and loosening connections.
2. Rinse equipment in warm water before placing into washer.
3. Load device into the washer and place small parts in container or tray inside the washer unit in order to avoid losing small components. **NOTE:** Load devices carefully into washer in order to avoid collision.
4. Follow the instructions listed below and set washer machine to these exact parameters:

Phase	Time (Min.)	Water Temperature	Detergent and Concentration
Pre-wash 1	4:00	Cold water	N/A
Enzyme Wash	4:00	Hot water	Neutral pH enzymatic (e.g. Endozime® AW Triple Plus with APA, Ratio 1:128)
Wash 1	10:00	60.0°C (140° F)	Neutral pH detergent (e.g. Renu-Klenz™, Steris, Ratio 1:256)
Rinse 1	0:30	Hot water	N/A
Thermal Rinse**	2:00	82.2°C (180° F)	N/A

\*\* Optional phase for disinfection of components – minimum water temperature as indicated or per worker manufacturer specifications for the thermal rinse cycle.

**NOTE:** Any deviation from this guideline could result in damage to the equipment as well as improper cleaning results. Rinse with purified water. Do not perform if parameters cannot be achieved.

5. Remove from washer and dry completely if needed.
6. Inspect equipment to make sure there is no visible organic debris or residue from cleaning agent. Repeat process if any visible soil is detected.

The life expectancy of the MAYFIELD products is expected to be as long as 7 years.

### **Care and Storage**

The MAYFIELD Adult Radiolucent Horseshoe Headrest should be returned to its protective storage case after each use. If the MAYFIELD Adult Radiolucent Horseshoe Headrest is dropped or mishandled, it should be immediately returned to Integra LifeSciences, Cincinnati, Ohio, for evaluation.

**WARNING:**

Device must not be used if it has been dropped or mishandled. Doing so may result in serious patient injury.

### **Maintenance and Care**

To ensure proper function and to extend the life and performance of the equipment, Integra LifeSciences recommends the following:

Recommended Action	Recommended Frequency
Return the device to the Integra LifeSciences Repairs department for detailed inspection and servicing	Once / year
Request that Integra NeuroSpecialists perform routine inspections of the device	Twice / year

In the absence of proper care and servicing of the device, negative effects may be seen after repeated processing over time which may lead to reduced performance.

Contact information: See the Service and Repair section for contact information on how to return your device for periodic servicing and to request periodic inspections.

See Inspection and/or Service notes section for routine checks to be performed on the device.

**NOTE:** Any serious incident that has occurred in relation to the device for the user and/or the patient should be reported to the manufacturer and the competent authority of the member state in which the user and/or patient is established.

### **Device Disposal**

**NOTE:** Follow hospital procedures for disposal of this device.

### **Integra Standard Warranty**

INTEGRA LIFESCIENCES CORPORATION (“INTEGRA”) warrants to the original purchaser only that each new MAYFIELD product is free from manufacturing defects in material and workmanship under normal use and service for a period of one year (except as otherwise expressly provided as to accessory items) from the date of delivery by INTEGRA to the first purchaser, but in no event beyond the expiration date stated on any product labeling.

- Surgical instruments are guaranteed to be free from defects in material and workmanship when maintained and cleaned properly and used normally for their intended purpose.
- Any covered product that is placed by INTEGRA under a lease, rental or installment purchase agreement and that requires repair service during the term of such placement agreement shall be repaired in accordance with the terms of such agreement.

If any covered defect occurs during the warranty period or term of such placement agreement, the purchaser should communicate directly with INTEGRA's home office. If purchaser seeks to invoke the terms of this warranty, the product must be returned to INTEGRA at its home office. The defective product should be returned promptly, properly packaged and postage prepaid. Loss or damage in return shipment to INTEGRA shall be at CUSTOMER's risk. INTEGRA's sole responsibility under this warranty shall be repair or replacement, at INTEGRA's sole discretion at INTEGRA's expense, subject to the terms of this warranty and applicable agreements.

IN NO EVENT SHALL INTEGRA BE LIABLE FOR ANY INCIDENTAL, INDIRECT, CONSEQUENTIAL OR PUNITIVE DAMAGES IN CONNECTION WITH THE ACQUISITION OR USE OF ANY INTEGRA PRODUCT. Further, this warranty shall not apply to, and INTEGRA shall not be responsible for, any loss arising in connection with the purchase or use of any INTEGRA product that has been repaired by anyone other than an authorized INTEGRA service representative or altered in any way so as, in INTEGRA's judgment, to affect its stability or reliability, or which has been subject to misuse, negligence or accident, or which has been used otherwise than in accordance with the instructions furnished by INTEGRA. THIS LIMITED WARRANTY IS EXCLUSIVE AND IN LIEU OF ALL OTHER WARRANTIES, EXPRESSED OR IMPLIED, AND OF ALL OTHER OBLIGATIONS OR LIABILITIES ON INTEGRA'S PART, AND INTEGRA NEITHER ASSUMES NOR AUTHORIZES ANY REPRESENTATIVE OR OTHER PERSON TO ASSUME FOR IT ANY OTHER LIABILITY IN CONNECTION WITH INTEGRA'S PRODUCTS.

INTEGRA DISCLAIMS ALL OTHER WARRANTIES, EXPRESSED OR IMPLIED INCLUDING ANY IMPLIED WARRANTY OF MERCHANTABILITY OR OF FITNESS FOR A PARTICULAR PURPOSE OR APPLICATION OR WARRANTY OF QUALITY AS WELL AS ANY EXPRESSED OR IMPLIED WARRANTY TO PATIENTS. No warranty or guarantee may be created by any act or statement nor may this Standard Warranty be modified in any way, except as a result of a writing signed by an officer of INTEGRA. These limitations on the creation or modification of this warranty may not be waived or modified orally or by any conduct.

### **Service and Repair**

For service and repairs outside the United States, contact your local authorized Integra representative.

Inside the United States, send all instruments for service or repair to:

Integra LifeSciences Corporation  
4900 Charlemer Drive, Building A  
Cincinnati, Ohio 45227

(Always include the purchase order number and a written description of the problem).

Or phone: 877-444-1114 (USA only)

Integra and the Integra logo are registered trademarks of Integra LifeSciences Corporation in the United States and/or other countries. MAYFIELD is a registered trademark of SM USA, Inc. and is used by Integra under license. Endozime is a trademark of Ruhof Corporation. Renu-Klenz is a trademark of Steris Corporation.

©2021 Integra LifeSciences Corporation. All Rights Reserved.  
451A2010 Rev DAA 2021-01 0687829-3

This page is intentionally left blank.

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