The AccuDrain® External CSF Drainage Systems

Sterile For Single Use Only

R_x ONLY

STERILE EO

CE 2797

INS-8400


INS-8401

Description

The AccuDrain® External CSF Drainage System is used to drain cerebrospinal fluid (CSF) from the ventricles of the brain or the lumbar subarachnoid space to a drainage bag. System components facilitate CSF drainage, CSF sampling and Intracranial Pressure (ICP) monitoring.

Monitoring of Intracranial Pressure (ICP) is usually performed in selected patients with severe head injury, subarachnoid hemorrhage, Reyes syndrome or similar encephalopathies, hydrocephalus, hydrocephalic shunt infections, intracranial hemorrhage and pre and/or post operative monitoring

Common features to the System include a closed sterile fluid path, calibrated 75ml graduated burette, 700ml drainage bag, anti-microbial hydrophobic vent, two needleless sampling sites, a patient stopcock, a four-way high flow burette stopcock, cord suspension, red burette height indicator, and integral pole attachment feature. A line level or laser device can be attached. The transparent graduated burette provides visualization of the flow and clarity of the CSF, and allows for volumetric readings to be made prior to fluid entering the drainage bag.

The AccuDrain External CSF Drainage Systems can be used with all of the Integra Ventricular and Lumbar Catheters. A tubing connector with a luer connector is provided to facilitate the catheter connection to the drainage system.

A fluid filled ICP monitoring transducer (not included) can be attached to the red capped port of the patient stopcock at the zero reference position.

Use of this device in an MR environment has been demonstrated to present no additional risk to the patient or other users.

The tubing, drainage bag and burette cap used in this product individually or in combination contains phthalates (DEHP). This device utilizes PVC tubing which is known to contain DEHP in a portion of the drainage path which on occasion is used for sampling and fluid injection. The risk of phthalates exposure to patients, including the male fetus, male neonate, and peri-pubescent male, is not considered significant for the following reasons: Sampled and drained CSF is not intended to be reintroduced into the patient. Injected fluids such as anti-coagulants and saline solution are only used as needed to eliminate or reduce blockage of the ventricular catheter and thus are not continuous in use.

Features and Benefits

The following features and benefits are common to each System.

- MRI Safe¹.
- This product does not contain and is not manufactured with dry natural rubber or natural rubber latex.
- Closed System
- Integrated Pressure Scale from

-5 to +30cm H₂O and -3 to +22mmHg, ± 1 scale unit.

- Anti-microbial Hydrophobic Vent
- Single-hand operated pole attachment feature
- Suspension cord locking feature for easy manipulation of System height
- Pole straps for additional security
- Pole Attachment Feature compatible with IV poles with diameters ranging from 3/4" to 1-1/8" (1.9 cm to 2.86 cm)
- Single-hand squeeze-lok™ tabs for easy & secure burette height adjustment
- Red burette height indicator to show the drainage pressure
- 75ml graduated burette in 1ml increments for accurate drainage monitoring and large bore outlet designed to resist clogging (1-19 ml: ± 1 ml, 20-75 ml: ± 2 ml)
- Removable 700ml drainage bag, with anti-reflux valve and drain tube, marked with approximate volume graduations from 100 to 700ml in 100ml ±20% increments
- Patient line comprised of green striped pressure tubing
- Line level, can be attached to zero reference
- Sterile packaged luer connector included for catheter connection
- Two needleless sampling sites

The AccuDrain External CSF Drainage System, INS-8400

This System includes all the features previously stated with the exception of an anti-reflux valve in the patient line.

The AccuDrain External CSF Drainage System, INS-8401

This System includes all the features previously stated with the addition of an integral low pressure anti-reflux valve in the patient line, distal to the zero reference stopcock.

Indications

Draining and monitoring of Cerebrospinal Fluid (CSF) flow from the ventricles of the brain or lumbar subarachnoid space is indicated in selected patients to:

- Reduce Intracranial Pressure (ICP)
- Monitor Intracranial Pressure (ICP)
- Monitor Cerebrospinal Fluid (CSF)
- Provide temporary CSF drainage

Contraindications

This device is not designed, sold, or intended for use except as indicated.

The External Drainage System is contraindicated in the following: Anticoagulation therapy, Coagulation disorders, Untreated scalp infections.

System use is contraindicated where trained personnel are not available to supervise drainage and monitoring on a 24-hour a day basis.

Warnings

¹The Laser level (INS-8902) is MR Unsafe. Removed it from the AccuDrain System prior to introducing the AccuDrain System into the MR environment.

Patients with cerebrospinal fluid drainage systems must be kept

under close observation for signs and symptoms of changing intracranial pressure. These signs and symptoms may vary from patient to patient. Increased intracranial pressure may be characterized by, but not limited to, headache, vomiting, irritability, listlessness, drowsiness, other signs of deterioration of consciousness and nuchal rigidity. In the infant, increased scalp tension at the anterior fontanelle and congestion of scalp veins may be noted.

Failure to appropriately adjust the rate of CSF outflow through the external drainage system may result in potentially serious injury to the patient. Improper drainage system setup can lead to overdrainage or underdrainage and potentially serious injury to the patient.

Proper alignment of the drainage system is critical whenever the system is open to drainage. When open, the height of the drainage system in relation to the patient affects the amount of CSF drainage and therefore, effects ICP. It is essential that neither the patient nor the drainage system be raised or lowered inadvertently. Whenever the device or patient is moved, the system should be realigned. Height changes should only be made by qualified personnel on the orders of the physician.

It is possible that the drainage system may lead to a false pressure reading. Reasons include (but are not limited to) the pressure line becoming clogged or kinked or from an air bubble lodged in the system. An incorrect pressure reading may lead to the wrong treatment for the patient.

In order to minimize the possibility of infection, meningitis or ventriculitis, the sampling site should be cleaned according to hospital protocol prior to use.

Strict aseptic technique should be used during initial set-up, at all stages of utilization and maintenance, and any time the system must be accessed, changed, or otherwise manipulated in any way that opens the closed system once it is connected to the patient.

If monitoring is not used, it is imperative to keep the red end cap located on the panel mount ICP monitoring port securely tightened.

Precautions

This product is intended for single use only. Reuse of the device can result in contamination and/or disease transmission.

This product should not be resterilized. Resterilization may affect the performance characteristics and the safety of the device.

Do not insert a needle into the needleless sampling site.

While the filter is hydrophobic, some manipulation may be required to reestablish drainage after repeated prolonged contact with CSF.

Integra makes no claim for or representation as to the performance characteristics of this product if it is used in conjunction with components of other manufacturers.

Contact with collodion remover, acetone and/or any other incompatible substances with product material may affect the integrity of the device, causing cracks and potential leaks on plastic components.

Complications

Complications which may result from the use of this product include the risks associated with the medi-

cations and methods utilized in the surgical procedure, as well as the patient's degree of intolerance to any foreign object implanted in the body.

The principal complications associated with cerebrospinal fluid drainage are infection, obstruction, or intracranial hypotension/hypertension.

The presence of a foreign body (i.e. the catheter system) may trigger infection or adverse reaction. Infection is a common and serious complication of a drainage system and is most frequently caused by skin contaminants. The incidence of these infections can be reduced by using aseptic technique. Risk of infection or other complications may be reduced by following hospital protocol for duration of use. Septicemia can result from infections anywhere in the body and may develop with few or no symptoms. Obstruction, partial obstruction, or irregular flow, may occur. Causes may include (but not limited to) obstruction by particulate matter (blood clots, fibrin, brain fragments), or malpositioning of the ventricular catheter.

Failure of the drainage system may be evidenced by continuing symptoms of increased ICP. Any failure of the system requires immediate replacement of the drainage system or the affected component.

Instructions for Use

Prior to system use, it is necessary for all responsible personnel to understand the use and function of the system components, system preparation, and system control.

Preparation and Alignment

The system should be prepared under sterile conditions prior to

placement of ventricular or lumbar catheter by a qualified health care professional.

Irrigation with preservative-free sterile normal saline may be performed prior to use as determined by hospital protocol. Caps may be temporarily loosened to allow air to escape. Check to ensure the absence of any residual air bubbles that may affect pressure transducer monitoring. Ensure that fluid flows from the burette to the drainage bag. When the panel mount ICP monitoring port is not in use, the red end cap must be in place to help prevent contamination.

1. Suspend the drainage system from the IV pole using the cord lock. Next, secure the drainage system to the IV pole by squeezing the pole attachment feature. Ensure that the The AccuDrain System is securely attached to the pole.
2. Once the System is secured, use the enclosed line level or laser level (not included) to ensure that the system has been properly adjusted with the zero reference on the pressure scale leveled to the external landmark as ordered by the physician, i.e. Foramen of Monro for ventricular catheters. Note: the line level or laser level is attached by inserting the fastener into the circular opening located on the back of the panel at the level of the zero reference point. Once attached, the system can be leveled with the line level by pulling the string taut to the anatomical landmark and aligning the bubble in the viewfinder. Once leveled, the line level can be attached by affixing the black o-ring to the hook located at the top back panel of the system. The system can be leveled with the laser level by aligning the bubble in the viewfinder and rotating the laser towards the

patient's anatomical landmark and then activating the laser to emit pointer beam. **Caution: Read the instructions for use (IFU) for the Integra laser level device before using laser.**

3. Tighten the suspension cord over the IV pole by moving the cord lock so that the suspension cord is taut. Secure the pole straps by pulling around the IV pole and inserting into the notch on each handle.
4. If an external transducer is to be used for monitoring ICP, calibration of the transducer should be completed, according to manufacturers instructions, prior to connecting to the patient. Hospital protocol should be used when using external transducers. **Caution: Method of transducer calibration may be affected by the presence of a one-way valve.**

To Set Pressure Level


Place the zero reference at the level of the external landmark ordered by the physician. Increase or decrease the height of the pressure level by manipulating the squeeze-lok to move the burette up or down relative to the patient.

Align the center of the red pressure indicator with the pressure level ordered (mmHg or cmH₂O) by the physician.

To Monitor Pressure

If closed system pressure monitoring and waveforms are desired, the AccuDrain System should be temporarily closed to drainage into the graduated burette. Turn stopcock OFF to burette to stop flow to the burette. **Note: reorient stopcock to resume drainage.**

To Monitor Pressure


To Drain CSF

Turn patient stopcock to allow CSF flow into the burette.

The ICP is affected by the volume of CSF drained and by the height of the pressure level of the burette. It is critical that neither the patient nor the drainage system be raised or lowered accidentally. Height changes should only be made by qualified personnel on the orders of a qualified physician.


To Drain CSF

To Collect and Measure CSF


Accurate determination of fluid accumulation may be accomplished with the graduated burette from 1ml to 75ml in 1ml increments. The burette stopcock should be OFF to the burette.

To Measure CSF Flow

To drain CSF from burette into the drainage bag, the stopcock should be OFF to the sampling site.

To Drain CSF from burette into the Drainage Bag

To Sample CSF

Sampling of CSF may be accomplished at the various sampling sites that have been cleaned according to hospital protocol. The SmartSite® needleless site can be accessed with any standard luer fitting syringe. Access may be accomplished at the patient line stopcock or the Y-site immediately distal to the patient stopcock or beneath the burette. If CSF access is desired from the burette, use the sampling site at the burette stopcock. **Do not insert a needle into the needleless sampling site.** In addition, CSF may be sampled from the bottom of the drainage bag using a syringe with a 25, or applicable gauge needle according to hospital protocol.

To Flush System

Sampling sites must be cleaned and flushed by qualified personnel, according to hospital protocol. Orient the stopcocks to temporarily prevent flow back to the patient. Ensure and verify that fluid flushes into the drainage bag.

Once flushing has been accomplished, reorient stopcocks to allow flow from the patient to the burette and into the drainage bag as ordered.

To Replace Drainage Bag

Orient the burette stopcock to temporarily stop flow into the drainage bag. Following hospital protocol clean the connection site between the burette and drainage bag.

Using aseptic technique, disconnect the drainage bag from the stopcock, attach tethered blue cap and detach bag from the system panel. Connect the sterile replacement drainage bag to drain tube. Attach bag to system panel.

Reorient burette stopcock to allow flow from the burette into the drainage bag. Verify flow into the drainage bag.

To Empty Drainage Bag

The drainage bag may be emptied using sterile technique, and a sterile syringe with a 25-gauge or applicable gauge needle.

To Transport Patient

If it is necessary to transport a patient while system is in full use, the system should remain vertical and correctly aligned with the patient for desired pressure level and drainage. If this is not possible, the graduated burette should be emptied into the drainage bag. The burette stopcock should then be

temporarily closed to prevent retrograde flow back into the burette. The panel mount stopcock should be oriented to temporarily stop flow from the patient to the burette.

CAUTION: If it becomes necessary to remove the system from the vertical secured position, orient the system burette side up and support the weight of the drainage bag when lifting or during movement.

The patient should then be transported as required and the system maintained according to physician orders and hospital protocol.

Note: During transport and procedures avoid excessive bending of the patient tubing line next to the three way stopcock assembly.

After patient transport has been completed, system use should be re-established with the correct pressure level reference relative to patient according to physician's orders. All stopcocks should be reoriented to reestablish flow. **Ensure and verify flow from the patient into the burette.**

DISPOSAL

After patient use, the systems must be handled as biohazardous material and disposed of in accordance with local and state environmental requirements following facility protocols.

Product Information Disclosure

INTEGRA HAS EXERCISED REASONABLE CARE IN THE CHOICE OF MATERIALS AND MANUFACTURE OF THIS PRODUCT. INTEGRA EXCLUDES ALL WARRANTIES, WHETHER EXPRESSED OR IMPLIED BY OPERATION OF LAW OR OTHERWISE, INCLUDING, BUT NOT LIMITED TO ANY IMPLIED WARRANTIES OF MERCHANTABILITY

OR FITNESS FOR A PARTICULAR PURPOSE. INTEGRA SHALL NOT BE LIABLE FOR ANY INCIDENTAL OR CONSEQUENTIAL LOSS, DAMAGE, OR EXPENSE, DIRECTLY OR INDIRECTLY ARISING FROM USE OF THIS PRODUCT. INTEGRA NEITHER ASSUMES OR AUTHORIZES ANY OTHER PERSON TO ASSUME FOR IT, ANY OTHER OR ADDITIONAL LIABILITY OR RESPONSIBILITY IN CONNECTION WITH THIS DEVICE.

Special Order Products

If this product is a Special Order product, i.e. requested by a physician, there may be differences between the enclosed product and the product description in this brochure. These differences will not affect the safety or effectiveness of the special order product.

How Supplied

Integra External CSF Drainage Systems are supplied sterile and non-pyrogenic in single wrap packaging.

Do Not Resterilize

All External CSF Drainage Systems are disposable devices. Integra does not recommend resterilization of these products.

Returned Goods Policy

Products must be returned in unopened packages, with manufacturer's seals intact to be accepted for replacement or credit, unless returned due to a complaint of product defect or mislabeling.

Determination of a product defect or mislabeling will be made by Integra, which determination will be final.

Products will be accepted for replacement if they have been in the possession of the customer for more than 90 days.

Storage Condition

Store at room temperature.

Product Order Information

All products can be ordered through your Integra Neuro Specialist or customer service representative or by contacting:

Integra NeuroSciences
1100 Campus Road
Princeton, NJ 08540 USA
Telephone: 1-800-654-2873
Outside the US: 609-275-0500
Fax: 609-275-5363
www.integralife.com

Integra LifeSciences Services (France)
SAS Immeuble Sequoia 2
97 allée Alexandre Borodine
Parc Technologique de la Porte des Alpes 69800 Saint Priest – France
Telephone: 33 (0) 4 37 47 59 10
Fax: 33 (0) 4 37 47 59 29

Caution: Federal (USA) law restricts this device to sale by or on the order of a physician. Do not use if the package has been opened or damaged.

Catalog Numbers

INS-8400 The AccuDrain External CSF Drainage System with Needleless sampling Sites, without anti-reflux valve

INS-8401 The AccuDrain External CSF Drainage System with Anti-Reflex Valve and Needleless sampling Sites

Accessory Catalog Numbers

INS-8700 Integra External CSF

Drainage System Replacement Bags /W Drain Tube. (MR Safe)
INS-8902 Integra Laser Level Device. (MR Unsafe)
INS-8903 Line Level Assembly (MR Safe)

Symbols Used On Labeling	
	Consult Instructions for Use
	Do not re-sterilize
	Do not use if package is damaged
	This product does not contain and is not manufactured with Dry Natural Rubber or Natural Rubber Latex
	Expiration date
	Do not re-use
	Lot number
	MR Safe
	Sterilized using Ethylene Oxide.
	Contains or Presence of DEHP
	Product complies with requirements of directive 93/42/EEC.
	Catalog number
	Manufacturer
	Authorized Representative in the European Community
Rx ONLY	Caution: Federal (USA) law restricts this device to sale by or on the order of a physician or practitioner.

Bibliography

1. Choi, S.C., Muizelaar, J.P. Barnes, T.Y., Maramarou, A., Brooks, D.M. and Young, H.F.: Prediction tree for severely head-injured patients. **J. Neurosurg.** 75:251-255; 1991.
2. Haan, J. and Thomeer, R.T.W.M.: Predictive value of temporary external lumbar drainage in normal pressure hydrocephalus. **Neurosurgery** 22: 388-391; 1988.
3. Ochiai H, Yamakawa Y., Continuous Lumbar Drainage for the Preoperative Management of Patients with Aneurysmal Subarachnoid Hemorrhage. **Neurol Med Chir (Tokyo)** 41: 576-581; 2001.
4. Roitberg, BZ, Khan, N, Alp, MS., Hersonskey, T., Charbell, FT., Ausman JI., Bedside External Ventricular Drain Placement for Treatment of Acute Hydrocephalus. **British Journal of Neurosurgery** : 15(4): 324-327; 2001
5. Scheinblum, S.T. and Hammond, M.: The treatment of Children with Shunt Infections: Extraventricular Drainage System Care. **Pediatric Nursing**, Vol. 16, 2: 139-143; 1990.
6. Thompson, HJ, Managing Patients with Lumbar Drainage Devices. **Critical Care Nurse**, Vol. 20, No.5; 2000.

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Bibliographie

- Choi, S.C., Muizelaar, J.P. Barnes, T.Y., Maramarou, A., Brooks, D.M. and Young, H.F.: Prediction tree for severely head-injured patients. *J. Neurosurg.* 75:251-255; 1991.
- Haan, J. and Thomeer, R.T.W.M.: Predictive value of temporary external lumbar drainage in normal pressure hydrocephalus. *Neurosurgery* 22: 388-391; 1988.
- Ochiai H, Yamakawa Y., Continuous Lumbar Drainage for the Preoperative Management of Patients with Aneurysmal Subarachnoid Hemorrhage. *Neurol Med Chir (Tokyo)* 41: 576-581; 2001.
- Roitberg, BZ, Khan, N, Alp, MS., Hersonskey, T., Charbell, FT., Ausman JI., Bedside External Ventricular Drain Placement for

Treatment of Acute Hydrocephalus. **British Journal of Neurosurgery** : 15(4): 324-327; 2001

5. Scheinblum, S.T. and Hammond, M.: The treatment of Children with Shunt Infections: Extraventricular Drainage System Care. **Pediatric Nursing**, Vol. 16, 2: 139-143; 1990.
6. Thompson, HJ, Managing Patients with Lumbar Drainage Devices. **Critical Care Nurse**, Vol. 20, No.5; 2000.

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Literatur

- Choi, S.C., Muizelaar, J.P. Barnes, T.Y., Maramarou, A., Brooks, D.M. and Young, H.F.: Prediction tree for severely head-injured patients. *J. Neurosurg.* 75:251-255; 1991.
- Haan, J. and Thomeer, R.T.W.M.: Predictive value of temporary external lumbar drainage in normal pressure hydrocephalus. *Neurosurgery* 22: 388-391; 1988.
- Ochiai H, Yamakawa Y., Continuous Lumbar Drainage for the Pre-operative Management of Patients with Aneurysmal Subarachnoid Hemorrhage. *Neurol Med Chir (Tokyo)* 41: 576-581; 2001.
- Roitberg, BZ, Khan, N, Alp, MS., Hersonskey, T., Charbell, FT., Ausman JI., Bedside External Ventricular Drain Placement for Treatment of Acute Hydrocephalus. *British Journal of Neurosurgery* : 15(4): 324-327; 2001
- Scheinblum, S.T. and Hammond, M.: The treatment of Children with Shunt Infections: Extraventricular

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Integra LifeSciences Services
(France)

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3. Ochiai H, Yamakawa Y., Continuous Lumbar Drainage for the Pre-operative Management of Patients with Aneurysmal Subarachnoid Hemorrhage. **Neurol Med Chir (Tokyo)** 41: 576-581; 2001.

4. Roitberg, BZ, Khan, N, Alp, MS., Hersonskey, T., Charbell, FT., Ausman JI., Bedside External Ventricular Drain Placement for Treatment of Acute Hydrocephalus. **British Journal of Neurosurgery** : 15(4): 324-327; 2001

5. Scheinblum, S.T. and Hammond, M.: The treatment of Children with Shunt Infections: Extraventricular Drainage System Care. **Pediatric Nursing**, Vol. 16, 2: 139-143; 1990.

6. Thompson, HJ, Managing Patients with Lumbar Drainage Devices. **Critical Care Nurse**, Vol. 20, No.5; 2000.

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- Choi, S.C., Muizelaar, J.P. Barnes, T.Y., Maramarou, A., Brooks, D.M. and Young, H.F.: Prediction tree for severely head-injured patients. *J. Neurosurg.* 75:251-255; 1991.
- Haan, J. and Thomeer, R.T.W.M.: Predictive value of temporary external lumbar drainage in normal pressure hydrocephalus. *Neurosurgery* 22: 388-391; 1988.
- Ochiai H, Yamakawa Y., Continuous Lumbar Drainage for the Pre-operative Management of Patients with Aneurysmal Subarachnoid Hemorrhage. *Neurol Med Chir (Tokyo)* 41: 576-581; 2001.
- Roitberg, BZ, Khan, N, Alp, MS., Hersonskey, T., Charbell, FT., Ausman JI., Bedside External Ventricular Drain Placement for Treatment of Acute Hydrocephalus. *British Journal of Neurosurgery* : 15(4): 324-327; 2001
- Scheinblum, S.T. and Hammond, M.: The treatment of Children with Shunt Infections: Extraventricular Drainage System Care. *Pediatric Nursing, Vol. 16*, 2: 139-143; 1990.
- Thompson, HJ, Managing

Patients with Lumbar Drainage
Devices. *Critical Care Nurse*, Vol.
20, No.5; 2000.

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Bibliografie

- Choi, S.C., Muizelaar, J.P. Barnes, T.Y., Maramarou, A., Brooks, D.M. and Young, H.F.: Prediction tree for severely head-injured patients. *J. Neurosurg.* 75:251-255; 1991.
- Haan, J. and Thomeer, R.T.W.M.: Predictive value of temporary external lumbar drainage in normal pressure hydrocephalus. *Neurosurgery* 22: 388-391; 1988.
- Ochiai H, Yamakawa Y., Continuous Lumbar Drainage for the Pre-operative Management of Patients with Aneurysmal Subarachoid Hemorrhage. *Neurol Med Chir (Tokyo)* 41: 576-581; 2001.
- Roitberg, BZ, Khan, N, Alp, MS., Hersonskey, T., Charbell, FT., Ausman JI., Bedside External Ventricular Drain Placement for Treatment of Acute Hydrocephalus. *British Journal of Neurosurgery* : 15(4): 324-327; 2001
- Scheinblum, S.T. and Hammond, M.: The treatment of Children with Shunt Infections: Extraventricular

Drainage System Care. **Pediatric Nursing, Vol. 16**, 2: 139-143; 1990.

6. Thompson, HJ, Managing Patients with Lumbar Drainage Devices. **Critical Care Nurse, Vol. 20**, No.5; 2000.

(This Page Left Intentionally Blank)

 Integra NeuroSciences
1100 Campus Road, Princeton, NJ 08540 USA
www.integralife.com

AccuDrain, Integra and the Integra logo are registered trademarks of Integra LifeSciences Corporation or its subsidiaries in the United States and/or other countries. All other trademarks and trade names are the property of their respective owners. © 2020 Integra LifeSciences Corporation. All rights reserved.